BILANCIA MULTIFUNZIONE

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BODY FAT SCALE

GENERAL SAFETY INSTRUCTIONS

Read all instructions and save it for future reference.

1. The appliance may be used by children of 8 years or older and by people with reduced physical, sensory or mental abilities or with lack of experience and the necessary knowledge, as long as they are supervised or have been instructed in the safe use of the appliance and have understood the dangers involved. Children must not play with the device. The cleaning and maintenance must not be carried out by unsupervised children.
2. Never leave children unsupervised with the packaging materials or the device. There is a risk of suffocation and loss of life.
3. **ATTENTION!** Do not stand on the scales with wet feet, or if the platform of the scales is wet — **RISK OF SLIPPING!**
4. Keep the platform dry, as it becomes very slippery when wet.
5. Do not use the appliance on wet or slippery ground.
6. Do not stand on the edge of the appliance.
7. Do not drop the appliance or jump on it; such actions could damage it.
8. Do not place any objects on step on platform when the appliance not in use.
9. Do not immerse the appliance in water or other liquids.
10. The batteries must not be exposed to sources of

heat such as fire, the sun, or similar sources.

11. Note on the batteries:

- A. The ingestion of batteries can be extremely dangerous. Keep the batteries and the scale out of the reach of children! In the event that a battery is ingested by a child, consult a doctor immediately.
- B. Replace spent batteries before they are completely dead.
- C. Always replace all of the batteries, with batteries of the same type.
- D. The batteries must not be recharged, taken apart, thrown in the fire or short-circuited.
- E. Electrolyte leaking from the batteries can damage the appliance. If you do not intend to use the appliance for a long period of time, remove the batteries from the battery compartment.
- F. If the batteries are leaking, wear gloves and dry the battery compartment well with a cloth.
- G. Batteries may contain toxins which can be harmful for health and the environment. Always dispose of batteries in accordance with applicable legal provisions. Do not dispose of batteries along with household waste.

12. When you remove the batteries at the end of their life or when the appliance is damaged, dispose of them in accordance with current national regulations as described below.

- In the event that fluid leaks from the batteries,

remove all of the batteries, taking care that the leaked fluid does not come into contact with skin or clothes. If the battery fluid does come into contact with skin or clothes, wash your skin with water immediately. Before inserting new batteries, clean the battery compartment with a damp cloth and allow it dry completely.

13. The batteries must be inserted correctly into the battery compartment, note the positive “+” and negative “-” sign. **WARNING:** Risk of explosion if the battery is replaced with an incorrect type. Use and replace only batteries of the same type.
14. Incorrect use of the batteries can cause fluid leakage, overheating or explosion. This fluid is corrosive and may be toxic. It may cause burns to skin and eyes, and is harmful if ingested.
15. To reduce the risk of accidents:
 - Keep the batteries out of the reach of children.
 - Do not heat up, open, pierce or cut the batteries or throw them into the fire.
 - Do not mix old and new batteries or batteries of different types.
 - Do not allow metal objects to come into contact with the terminals of the battery, as they may overheat and/or cause burns.


INTENDED USE

This product is intended exclusively for personal use indoors, and is not suitable for commercial or medical purposes. Use this scale you will be able to measure below data at home:

- Weight
- Body fat
- Water percentage
- Muscle mass
- Bone mass
- Calories
- Body mass index

NOTE: if need more precise data, it is better to check with doctor.

PACKAGE CONTENTS AND DESCRIPTION OF THE APPLIANCE

1. Large LCD screen
2. Select UP touch button
3. SET touch button
4. Select DOWN touch button
5. Step on platform
6. 2 high precision sensor bars
7. Battery compartment
8. Unit button (kg or lb)
9. Anti-slip feet


BEFORE THE FIRST USE

1. Remove all packaging materials and any protective films from the appliance.
2. Clean the step on platform with damp cloth.
3. Allow the step on platform dry completely before using it.

NOTE: the MAX weight measurement up to 180kg.

HOW TO USE

INSERT THE BATTERIES AND SET THE UNIT OF MEASUREMENT (KG/LB)

1. Remove the battery compartment's cover.
2. Follow the guidance of battery compartment, insert 2 AAA batteries, note the positive "+" and negative "-" signs on battery.
3. Place the cover back on to the battery compartment securely.
4. Once the batteries are inserted, press the unit button at backside of the scale to set the unit of measurement of the scale (kg or lb).

NOTE: when **LO** appear on the display screen it indicating the batteries needs to be replaced.

CALIBRATE THE SCALE AND SIMPLY WEIGHING

1. In order to use the scale it must be calibrated first, place the scale on a flat surface on the floor.

2. Step on the platform, the scale will turn on.
3. Continue to stand on the scale for a few moments while the scale reads your weight.
4. Your weight will appear on the display screen.
5. Once you step off the platform, the scale will auto turn off.

NOTE: the MAX weight measurement up to 180kg.

SET USER DATA


- Create a new user and enter user's personal information into scale.
- The scale can store up to 10 users.

1. Activate the scale, when the "0.0" appears on display screen, and then touch the  button. 

2. A "0" will start flash on the left side of screen indicating the number of user. 

3. **SET USER NUMBER:** touch the   buttons to set a number between 0-9, touch the  button again to save it. 

4. **SET GENDER:** once the user number is set, the gender icon start flashing, touch the   buttons to set the gender, and then touch the  button again to save it. 

5. **SET AGE:** once the user gender is set, the Age start flashing, touch the   buttons to set user's age, and then touch the  button again to save it. 

6. **SET HEIGHT:** once the user age is set, the Height start flashing, touch the   buttons to set user's Height, and then touch the  button again to save it. 


7. When the user's personal data is set. User number, user gender and "0.0" will appear on display screen indicating the scale is ready for advanced measurement. 

NOTE: if you make a mistake before finishing the data input, repeat the operations described above starting from the selection of the personal data code.


ADVANCED MEASUREMENT

- To guarantee the accuracy of the data reading, step onto the scales in your bare feet.
- Ensure that your feet are correctly aligned with the 2 high precision sensor bars on the platform when weighing.

1. Activate the scale, when the "0.0" appears on display screen, and then touch the  button once, the user number will start flash, touch the   buttons to select user number who is going to weighing, and the user gender will following to appear on left side of screen.

2. When the user number, user gender and "0.0" appear steady, the scale is ready for advanced measurement. 

3. Your weight value will be displayed and flash 3 times after stepping on the scale. Then, you will see "0" appearing on the scale, and then body fat, water percentage, muscle mass, bone mass, calories, body mass index will be displayed.

4. Scale will auto turn off after the measurement data displayed three times.

BODY MASS INDEX CHART (BODY FAT PERCENTAGE)

Age	Female			
	Body fat			
	Low	Healthy	High	Very high
12-20	<18%	18-28%	28-33%	>33%
21-42	<20%	20-30%	30-35%	>35%
43-65	<21%	21-31%	31-36%	>36%
66-100	<22%	22-32%	32-37%	>37%
Age	Male			
	Body fat			
	Low	Healthy	High	Very high
12-20	<15%	15-21%	21-26%	>26%
21-42	<17%	17-23%	23-28%	>28%
43-65	<18%	18-24%	24-29%	>29%
66-100	<19%	19-25%	25-30%	>30%

BODY MASS INDEX CHART (WATER HYDRATION)

Age	Male	Female	Index condition	Capacity: 37.8-66.0% Water percentage division: 0.1%
<=30	66.0-60.5%	66.0-59.9%	Very high	
	60.4-57.1%	59.8-56.4%	High	
	57.0-53.6%	56.3-53.0%	Healthy	
	53.5-50.2%	52.9-48.6%	Low	
	50.1-37.8%	48.5-37.8%	Very low	
>30	66.0-59.1%	66.0-55.0%	Very high	
	59.0-55.7%	54.9-51.6%	High	
	55.6-52.3%	51.5-48.1%	Healthy	
	52.2-48.8%	48.0-44.7%	Low	
	48.7-37.8%	44.6-37.8%	Very low	

BODY MASS INDEX CHART (MUSCLE MASS PERCENTAGE)

Age	Female			
	Muscle mass			
	Low	Healthy	High	Very high
18-40	<24.5	24.5-30.5%	30.6-35.5%	>35.5%
41-60	<24.3	24.3-30.3%	30.4-35.3%	>35.3%
61-80	<24.1	24.1-30.1%	30.2-35.1%	>35.1%
Age	Male			
	Muscle mass			
	Low	Healthy	High	Very high
18-40	<33.5	33.5-39.5%	39.6-45.0%	>45%
41-60	<33.3	33.3-39.3%	39.4-44.8%	>44.8%
61-80	<33.1	33.1-39.1%	39.2-44.6%	>44.6%

BODY MASS INDEX CHART (WEIGHT AND BONE MASS IN KILOGRAM)

Female	
Weight	Bone mass
0-52kg	2kg
52-77kg	2.5kg
77kg and up	3kg
Male	
Weight	Bone mass
0-63kg	2.7kg
63-97.5kg	3.4kg
97.5kg and up	3.6kg

BODY MASS INDEX CHART (CALORIES)

Age	Calories	
	Female	Male
10 ~ 40	655 ~ 3999	169 ~ 3603
41 ~ 56	655 ~ 3949	66 ~ 3392
57 ~ 85	655 ~ 3874	66 ~ 3283

BODY MASS INDEX CHART

Under weight	Healthy weight	Overweight
<18.5	18.5-25	>25

TECHNICAL DATA

1. Product size: 325x325mm
2. Max. Weight: 180kg
3. Unit: kg/lb
4. Over load reminder
5. Low battery reminder
6. Battery: AAA x2
7. Sensor: 2 high precision sensor bars
8. User: up to 10 users

CLEAN AND CARE

- Do not immerse any part of the appliance in water or other liquids.
- Clean the step on platform gently with a soft, damp cloth.
- Do not use aggressive, caustic, abrasive, and corrosive or brushes with stiff bristles for cleaning.
- Protect the appliance from damp, dust, chemicals, impacts and sources of heat.
- Do not place any objects on top of the scale when it is not in use.
- Store the appliance in a cool and dry place.

WARRANTY

The appliance is guaranteed for two years from the date of delivery, as shown by the date on the receipt/invoice. In case of defect existing before purchase, the replacement is guaranteed. The warranty does not cover those parts found to be faulty as a result of negligence or careless use. The warranty also excludes all cases of improper use of the appliance or any instances of use of professional type. The seller refuses to accept any liability for damage that may directly or indirectly be caused to people, property or pets as a result of failure to abide by all instructions contained in this manual.

CORRECT DISPOSAL OF THIS PRODUCT

The symbol representing the crossed out trash can on the appliance indicates that the product, at the end of its useful life, must be treated separately from domestic waste. It must be taken to a recycling centre  for electrical and electronic appliances or turned in to the point of sale when purchasing a new, equivalent, appliance. Correct consignment of the no-longer used appliance to a recycling centre for treatment and environmentally compatible elimination avoids negative impact on the environment and public health, and makes possible the recovery the materials that the product is made from. The user is responsible for taking the appliance to the recycling centre at the end of its life. For more detailed information regarding the available collection systems, contract the local refuse centre or the store where the product was purchased. Anyone who leaves or throws away this appliance without taking it to a collection centre for electrical-electronic waste, shall be punished by administrative fine as per applicable law on illegal waste disposal.

Exhausted batteries must also be disposed of in compliance with the environmental regulations in force, at an authorized collection center, or deposited in the appropriate containers at the point of sale where the purchase was made.

Batteries must be removed from the appliance before its disposal.

Do not throw spent or damaged batteries into household waste so as not to harm the environment.

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SPINSERVICE S.r.l.